

JOBST® FarrowWrap®

A-D Below Knee

Ready-to-wear order form and prescription request

Fax order to Customer Services on: **0845 122 3450**

Email order to Customer Services on: **compression.uk@jobst.com**

Date: _____ Purchase Order No.: _____

Patient Name: _____ DoB: _____

Measured By: _____

Tel: _____ Email: _____

Delivery Address: _____

Invoice Address: _____





Measuring instructions

Legpiece - 3 simple measuring points

1. Measure the widest calf **C** and narrowest ankle **B** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column (if 2 measurements do not fit in same column the smaller size should accommodate reducing the oedema)
3. Measure the length of the leg posteriorly from the floor to **D** measurement (2cm below the crease of the knee). Follow the contour of the leg. Use the sizing chart to select regular or tall to determine size

Footpiece - 2 simple measuring points

1. Measure **A'** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column
3. Measure the **X** length of the foot. Use the sizing chart to select regular or long to determine size



		X Small	Small	Medium	Large	X Large
Legpiece						
C	Mid-Calf	36-43cm	42-50cm	48-58cm	53-63cm	58-68cm
B	Ankle	21-25cm	25-30cm	30-36cm	36-42cm	42-50cm
A-D	Regular	35-39cm	37-41cm	39-43cm	41-45cm	42-46cm
A-D	Tall	40-44cm	42-46cm	44-48cm	46-50cm	47-51cm
NOTE: for A-D measure posteriorly and follow contours						
Footpiece						
A'	Mid-foot	22-24cm	25-27cm	28-30cm	31-34cm	35-40cm
X	Regular	15-17cm	16-19cm	18-20cm	19-22cm	21-24cm
X	Tall	18-20cm	20-22cm	21-23cm	23-25cm	25-27cm
NOTE: for X measure from the heel to the base of the little toe						

Note: 3cm allowance on either side of measurement.

Prescription request

If requesting a prescription please complete order form and return to the prescriber. The codes for specific products can be found on this form

Legpiece

Size: X Small Small Medium Large X Large

Length A-D: Regular Tall

Fabric: Lite Classic Strong

PIP Code: -

Code: -

Quantity Required:

Repeat prescription required every months

Footpiece

Size: X Small Small Medium Large X Large

Length X: Regular Long

Fabric: Lite Classic Strong

PIP Code: -

Code: -

Quantity Required:

Repeat prescription required every months

JOBST FarrowHybrid

Size: Small Medium Large X Large

Width: Regular* Wide**

Compression: 20-30mmHg - Taupe Black*** 30-40mmHg - Sand

PIP Code: -

Code: -

Quantity Required:

Repeat prescription required every months

*Regular: 20-27 cm **A'** **Wide: 28-35 cm **A'** ***Available in wide variant only

A pair of appropriately sized liners will be provided free of charge with each JOBST FarrowWrap lower limb piece.

Customer Service: **0845 122 3600** Compression Therapy Helpline: **0800 389 8424** Email: **compression.uk@jobst.com** Website: **www.bsnmedical.co.uk**

JOBST FarrowWrap Order Codes

JOBST FarrowWrap Lite (20-30mmHg)

Code	Size	Length	Colour	PIP Code
Legpiece				
76654-00	Extra Small	Regular	Beige	360-6233
76654-05	Extra Small	Tall	Beige	360-6282
76654-01	Small	Regular	Beige	360-6241
76654-06	Small	Tall	Beige	360-6290
76654-02	Medium	Regular	Beige	360-6258
76654-07	Medium	Tall	Beige	360-6308
76654-03	Large	Regular	Beige	360-6266
76654-08	Large	Tall	Beige	360-6316
76654-04	Extra Large	Regular	Beige	360-6274
76654-09	Extra Large	Tall	Beige	360-6324
Footpiece				
76652-05	Extra Small	Regular	Beige	360-6332
76652-00	Extra Small	Long	Beige	360-6381
76652-06	Small	Regular	Beige	360-6340
76652-01	Small	Long	Beige	360-6399
76652-07	Medium	Regular	Beige	360-6357
76652-02	Medium	Long	Beige	360-6407
76652-08	Large	Regular	Beige	360-6365
76652-03	Large	Long	Beige	360-6415
76652-09	Extra Large	Regular	Beige	360-6373
76652-04	Extra Large	Long	Beige	360-6423



JOBST FarrowWrap Classic (30-40mmHg)

Code	Size	Length	Colour	PIP Code
Legpiece				
76657-00	Extra Small	Regular	Beige	360-6589
76657-05	Extra Small	Tall	Beige	360-6639
76657-01	Small	Regular	Beige	360-6597
76657-06	Small	Tall	Beige	360-6647
76657-02	Medium	Regular	Beige	360-6605
76657-07	Medium	Tall	Beige	360-6654
76657-03	Large	Regular	Beige	360-6613
76657-08	Large	Tall	Beige	360-6662
76657-04	Extra Large	Regular	Beige	360-6621
76657-09	Extra Large	Tall	Beige	360-6670
Footpiece				
76656-05	Extra Small	Regular	Beige	360-6688
76656-00	Extra Small	Long	Beige	360-6738
76656-06	Small	Regular	Beige	360-6696
76656-01	Small	Long	Beige	360-6746
76656-07	Medium	Regular	Beige	360-6704
76656-02	Medium	Long	Beige	360-6753
76656-08	Large	Regular	Beige	360-6712
76656-03	Large	Long	Beige	360-6761
76656-09	Extra Large	Regular	Beige	360-6720
76656-04	Extra Large	Long	Beige	360-6779



JOBST FarrowWrap Strong (30-40mmHg)

Code	Size	Length	Colour	PIP Code
Legpiece				
76648-00	Extra Small	Regular	Beige	382-0073
76648-05	Extra Small	Tall	Beige	382-0024
76648-01	Small	Regular	Beige	382-0107
76648-06	Small	Tall	Beige	382-0057
76648-02	Medium	Regular	Beige	382-0099
76648-07	Medium	Tall	Beige	382-0040
76648-03	Large	Regular	Beige	382-0081
76648-08	Large	Tall	Beige	382-0032
76648-04	Extra Large	Regular	Beige	382-0065
76648-09	Extra Large	Tall	Beige	382-0016
Footpiece				
76646-04	Extra Small	Regular	Beige	381-9976
76646-09	Extra Small	Long	Beige	381-9927
76646-05	Small	Regular	Beige	382-0008
76646-00	Small	Long	Beige	381-9950
76646-06	Medium	Regular	Beige	381-9992
76646-01	Medium	Long	Beige	381-9943
76646-07	Large	Regular	Beige	381-9984
76646-02	Large	Long	Beige	381-9935
76646-08	Extra Large	Regular	Beige	381-9968
76646-03	Extra Large	Long	Beige	381-9919

JOBST FarrowHybrid* - supplying compression to the foot and ankle

Code	Size	Colour	PIP Code	Code	Size	Colour	PIP Code
ADI (20-30mmHg)				ADII (30-40mmHg)			
76664-00	Small (shoe size 2-5)	Beige	408-1626	76664-04	Small (shoe size 2-5)	Beige	408-1691
76664-01	Medium (shoe size 6-8)	Beige	408-1634	76664-05	Medium (shoe size 6-8)	Beige	408-1709
76664-02	Large (shoe size 9-11)	Beige	408-1642	76664-06	Large (shoe size 9-11)	Beige	408-1717
76664-03	Extra Large (shoe size 12-13)	Beige	408-1659	76664-07	Extra Large (shoe size 12-13)	Beige	408-1725
76664-12	Small - wide (shoe size 2-5)	Beige	408-1667	76664-08	Small - wide (shoe size 2-5)	Beige	408-1733
76664-13	Medium - wide (shoe size 6-8)	Beige	408-1675	76664-09	Medium - wide (shoe size 6-8)	Beige	408-1741
76664-14	Large - wide (shoe size 9-11)	Beige	408-2384	76664-10	Large - wide (shoe size 9-11)	Beige	408-1758
76664-15	Extra Large - wide (shoe size 12-13)	Beige	408-1683	76664-11	Extra Large - wide (shoe size 12-13)	Beige	408-1766
76664-16	Small - wide (shoe size 2-5)	Black	408-6765	*Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm. Consider a JOBST FarrowWrap footpiece in this circumstance.			
76664-17	Medium - wide (shoe size 6-8)	Black	408-6757				
76664-18	Large - wide (shoe size 9-11)	Black	408-6740				
76664-19	Extra Large - wide (shoe size 12-13)	Black	408-6732				



By completing this order form, you are confirming that you are aware of your obligation to obtain informed written consent from the patient on the processing of their data for the production of their JOBST compression garment in accordance with the General Data Protection Regulation. For more information and our Privacy Policy, visit www.bsnmedical.co.uk