

JOBST® FarrowWrap®

A-G Thigh High

Custom-fit / Direct purchasing order form


Fax order to customer services on: **0845 122 3450**

Email order to customer services on: **compression.uk@bsnmedical.com**

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____



Please measure all points for a JOBST® FarrowWrap® custom-fit product



E-G

Left: _____

Right: _____

D-E

Left: _____

Right: _____

Required for thighpiece

A-D

Left: _____

Right: _____

Note:
For A-D and D-E, measure posteriorly and follow leg contour.

Note:
For E-G measure medially and follow contours.

Custom-fit

Thighpiece A-G Lite Classic Strong

Legpiece A-D Lite Classic Strong

Footpiece A¹-X Lite Classic Strong

Quantity Required:

Repeat prescription required every months

Remarks

Please measure all points below for a JOBST® FarrowWrap® ready-to-wear product

	X Small	Small	Medium	Large	X Large	XX Large
Thighpiece (including kneepiece)						
G Groin	60-70 cm	65-75 cm	70-80 cm	75-90 cm	85-100 cm	-
E Above knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm	-
E-G Short	16-20 cm	16-20 cm	16-20 cm	16-20 cm	16-20 cm	-
E-G Regular	21-26 cm	21-26 cm	21-26 cm	21-26 cm	21-26 cm	-
E-G Tall	27-31 cm	27-31 cm	27-31 cm	27-31 cm	27-31 cm	-
NOTE: for E-G measure medially and follow contours						
Kneepiece (seperate)						
E Above knee	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm	81-87 cm
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm	63-73 cm
Legpiece						
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm	-
B Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm	-
A-D Regular	35-39 cm	37-41 cm	39-43 cm	41-45 cm	42-46 cm	-
A-D Tall	40-44 cm	42-46 cm	44-48 cm	46-50 cm	47-51 cm	-
NOTE: for A-D measure posteriorly and follow contours						
Footpiece						
A¹ Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm	-
X Regular	15-17 cm	16-19 cm	18-20 cm	19-22 cm	21-24 cm	-
X Tall	18-20 cm	20-22 cm	21-23 cm	23-25 cm	25-27 cm	-
NOTE: for X measure from the heel to the base of the little toe						

Note: 3cm allowance on either side of measurement.

Please use ready-to-wear product code required found on the back of this form

Ready-to-wear

Thighpiece -

Legpiece -

Footpiece -

Quantity Required:

Repeat prescription required every months

JOBST FarrowHybrid
 -

A-D Below Knee liner
 -

A-G Thigh High liner
 -

Remarks

Note: Custom-fit is not available on FP10 / GP10. A pair of appropriately sized liners will be provided free of charge with each JOBST FarrowWrap lower limb piece.

JOBST FarrowWrap ready-to-wear direct order form

JOBST FarrowWrap Lite (20-30mmHg)

Code	Size	Length
Thighpiece (including kneepiece)		
76655-05	Extra Small	Short
76655-00	Extra Small	Regular
76655-10	Extra Small	Tall
76655-06	Small	Short
76655-01	Small	Regular
76655-11	Small	Tall
76655-07	Medium	Short
76655-02	Medium	Regular
76655-12	Medium	Tall
76655-08	Large	Short
76655-03	Large	Regular
76655-13	Large	Tall
76655-09	Extra Large	Short
76655-04	Extra Large	Regular
76655-14	Extra Large	Tall
Kneepiece (seperate)*		
73367-00	Extra Small	-
73367-01	Small	-
73367-02	Medium	-
73367-03	Large	-
73367-04	Extra Large	-
Legpiece		
76654-00	Extra Small	Regular
76654-05	Extra Small	Tall
76654-01	Small	Regular
76654-06	Small	Tall
76654-02	Medium	Regular
76654-07	Medium	Tall
76654-03	Large	Regular
76654-08	Large	Tall
76654-04	Extra Large	Regular
76654-09	Extra Large	Tall
Footpiece		
76652-05	Extra Small	Regular
76652-00	Extra Small	Long
76652-06	Small	Regular
76652-01	Small	Long
76652-07	Medium	Regular
76652-02	Medium	Long
76652-08	Large	Regular
76652-03	Large	Long
76652-09	Extra Large	Regular
76652-04	Extra Large	Long



JOBST FarrowWrap Strong (30-40mmHg)

Code	Size	Length
Thighpiece (including kneepiece)		
76649-05	Extra Small	Short
76649-00	Extra Small	Regular
76649-10	Extra Small	Tall
76649-06	Small	Short
76649-01	Small	Regular
76649-11	Small	Tall
76649-07	Medium	Short
76649-02	Medium	Regular
76649-12	Medium	Tall
76649-08	Large	Short
76649-03	Large	Regular
76649-13	Large	Tall
76649-09	Extra Large	Short
76649-04	Extra Large	Regular
76649-14	Extra Large	Tall
Kneepiece (seperate)*		
73368-00	Extra Small	-
73368-01	Small	-
73368-02	Medium	-
73368-03	Large	-
73368-04	Extra Large	-
73368-05	Extra Extra Large	-
Legpiece		
76648-00	Extra Small	Regular
76648-05	Extra Small	Tall
76648-01	Small	Regular
76648-06	Small	Tall
76648-02	Medium	Regular
76648-07	Medium	Tall
76648-03	Large	Regular
76648-08	Large	Tall
76648-04	Extra Large	Regular
76648-09	Extra Large	Tall
Footpiece		
76646-04	Extra Small	Regular
76646-09	Extra Small	Long
76646-05	Small	Regular
76646-00	Small	Long
76646-06	Medium	Regular
76646-01	Medium	Long
76646-07	Large	Regular
76646-02	Large	Long
76646-08	Extra Large	Regular
76646-03	Extra Large	Long



JOBST FarrowWrap Classic (30-40mmHg)


Code	Size	Length
Thighpiece (including kneepiece)		
76658-06	Extra Small	Short
76658-00	Extra Small	Regular
76658-12	Extra Small	Tall
76658-07	Small	Short
76658-01	Small	Regular
76658-13	Small	Tall
76658-08	Medium	Short
76658-02	Medium	Regular
76658-14	Medium	Tall
76658-09	Large	Short
76658-03	Large	Regular
76658-15	Large	Tall
76658-10	Extra Large	Short
76658-04	Extra Large	Regular
76658-16	Extra Large	Tall
Legpiece		
76657-00	Extra Small	Regular
76657-05	Extra Small	Tall
76657-01	Small	Regular
76657-06	Small	Tall
76657-02	Medium	Regular
76657-07	Medium	Tall
76657-03	Large	Regular
76657-08	Large	Tall
76657-04	Extra Large	Regular
76657-09	Extra Large	Tall
Footpiece		
76656-05	Extra Small	Regular
76656-00	Extra Small	Long
76656-06	Small	Regular
76656-01	Small	Long
76656-07	Medium	Regular
76656-02	Medium	Long
76656-08	Large	Regular
76656-03	Large	Long
76656-09	Extra Large	Regular
76656-04	Extra Large	Long

Garment Grip Roll*

Code	Size
Garment Grip Roll	
76672-00	14cm x 80cm
76672-01	14cm x 120cm
76672-02	14cm x 550cm
76672-03	20cm x 120cm
76672-04	30cm x 150cm
76672-05	10cm x 60cm

JOBST FarrowHybrid[†] - supplying compression to the foot and ankle

Code	Size	Colour
ADI (20-30mmHg)		
76664-00	Small (shoe size 2-5)	Beige
76664-01	Medium (shoe size 6-8)	Beige
76664-02	Large (shoe size 9-11)	Beige
76664-03	Extra large (shoe size 12-13)	Beige
76664-12	Small - wide (shoe size 2-5)	Beige
76664-13	Medium - wide (shoe size 6-8)	Beige
76664-14	Large - wide (shoe size 9-11)	Beige
76664-15	Extra large - wide (shoe size 12-13)	Beige
76664-16	Small - wide (shoe size 2-5)	Black
76664-17	Medium - wide (shoe size 6-8)	Black
76664-18	Large - wide (shoe size 9-11)	Black
76664-19	Extra large - wide (shoe size 12-13)	Black
ADII (30-40mmHg)		
76664-04	Small (shoe size 2-5)	Beige
76664-05	Medium (shoe size 6-8)	Beige
76664-06	Large (shoe size 9-11)	Beige
76664-07	Extra large (shoe size 12-13)	Beige
76664-08	Small - wide (shoe size 2-5)	Beige
76664-09	Medium - wide (shoe size 6-8)	Beige
76664-10	Large - wide (shoe size 9-11)	Beige
76664-11	Extra large - wide (shoe size 12-13)	Beige



[†]Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm. Consider a JOBST FarrowWrap footpiece in this circumstance.

A-D Below Knee Liners*

Code	Size
Silver	
76668-00	Small (15-58cm)
76668-02	Large (52-70cm)
TG Soft	
76669-00	Small (<40cm)
76669-02	Medium (40-70cm)
76669-04	Large (70-120cm)

A-G High High Liners*

Code	Size
TG Soft	
76670-02	Medium (40-70cm)
76670-04	Large (70-120cm)

Foam Squares Pack*

Code	Size
Foam Squares Pack	
76675-01	Small (8x8cm)
76675-02	Large (10x10cm)

NB The thighpiece comes with kneepiece included, however, the kneepiece can be ordered separately.

*Not available on FP10 / GP10.

By completing this order form, you are confirming that you are aware of your obligation to obtain informed written consent from the patient on the processing of their data for the production of their JOBST compression garment in accordance with the General Data Protection Regulation. For more information and our Privacy Policy, visit www.bsnmedical.co.uk