


# JOBST® UlcerCare

## 2-in-1 Compression System

Ready-to-wear order form and prescription request

Email order to Customer Services on: **compression.uk@jobst.com**  
 Fax order to Customer Services on: **0345 122 3450**

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_  
 \_\_\_\_\_



### Measuring instructions

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Determine size (**S-4XL**) using chart below

Size	S	M	L	XL	2XL	3XL	4XL
<b>C</b>	29 - 35cm	34 - 41cm	40 - 47cm	46 - 54cm	48 - 56cm	50 - 58cm	52 - 60cm
<b>B</b>	18 - 21cm	21 - 24cm	24 - 27cm	27 - 30cm	30 - 33cm	33 - 36cm	36 - 39cm

**Top Tip**  
 If the measurements fall into two sizes (e.g. B is Medium and C is Large), choose the size with the measurement that is closest to the middle of the range.  
 For example, B is 22cm and C is 42cm.  
 22cm is in the middle of the Medium size (21-24cm).  
 42cm is at the start of the Large size (40-47cm).  
 In this example, we would advise you order the Medium garment.

### Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST® UlcerCare



**JOBST® UlcerCare Kit**

Each kit contains one outer stocking and two liners

**Size**

Small       Medium       Large  
 X Large       2X Large       3X Large  
 4X Large

**Colour**

Beige       Black

**Zipper** (to be placed on the side opposite to the ulceration)

No zip       Zip left       Zip right

PIP Code:    -

Product Code:      -

Quantity Required:

Repeat prescription required every  months

**JOBST® UlcerCare Liners**

Each pack contains three liners

**Size**

Small       Medium       Large  
 X Large       2X Large       3X Large  
 4X Large

PIP Code:    -

Product Code:      -

Quantity Required:

Repeat prescription required every  months

# JOBST® UlcerCare Order Codes

## JOBST® UlcerCare Kit (40mmHg)

- 1 x Outer Stocking
- 2 x Liners
- With Zip

Product Code	Size	PIP Code
<b>Zipper to be placed on the side opposite to the ulceration</b>		
<b>Zipper Left - Black</b>		
73631-41	S	364-1263
73631-42	M	364-1248
73631-43	L	364-1222
73631-44	XL	364-1206
73631-45	2XL	364-1172
73631-46	3XL	364-1156
73631-47	4XL	364-1131
<b>Zipper Right - Black</b>		
73631-71	S	364-1404
73631-72	M	364-1388
73631-73	L	364-1362
73631-74	XL	364-1347
73631-75	2XL	364-1321
73631-76	3XL	364-1305
73631-77	4XL	364-1289
<b>Zipper Left - Beige</b>		
73631-21	S	364-1271
73631-22	M	364-1255
73631-23	L	364-1230
73631-24	XL	364-1214
73631-25	2XL	364-1180
73631-26	3XL	364-1164
73631-27	4XL	364-1149
<b>Zipper Right - Beige</b>		
73631-61	S	364-1412
73631-62	M	364-1396
73631-63	L	364-1370
73631-64	XL	364-1354
73631-65	2XL	364-1339
73631-66	3XL	364-1313
73631-67	4XL	364-1297

## JOBST® UlcerCare Kit (40mmHg)

- 1 x Outer Stocking
- 2 x Liners
- Without Zip

Product Code	Size	PIP Code
<b>Black</b>		
73630-41	S	364-4556
73630-42	M	364-4549
73630-43	L	364-4531
73630-44	XL	364-4523
73630-45	2XL	364-4515
73630-46	3XL	364-4507
73630-47	4XL	364-4499
<b>Beige</b>		
73630-21	S	364-4622
73630-22	M	364-4614
73630-23	L	364-4606
73630-24	XL	364-4598
73630-25	2XL	364-4580
73630-26	3XL	364-4572
73630-27	4XL	364-4564

## JOBST® UlcerCare Liners (15-20mmHg)

- 3 x Liners

Product Code	Size	PIP Code
<b>Replacement Liners - White</b>		
73632-21	S	364-7666
73632-22	M	364-7674
73632-23	L	364-7682
73632-24	XL	364-7708
73632-25	2XL	364-1123
73632-26	3XL	364-1115
73632-27	4XL	364-1107

**Find ordering and prescribing information fast with JOBST Code Finder**

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for [codefinder.jobst.com](http://codefinder.jobst.com) on your mobile device.



**Fast and free home delivery of prescription garments with JOBST Delivered\***

To use this dispensing service, send this form and prescription to: **JOBST Delivered, FREEPOST, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG** or order online at [order.jobst.com/UK](http://order.jobst.com/UK) and select JOBST Delivered as the dispenser (note: prescription to be posted to JOBST Delivered unless using EPS).

If the patient is using an Electronic Prescription Service (EPS), select code **FJ708** to nominate Daylong Direct as their dispenser of choice.

For envelopes, please contact: [jobstdelivered@daylongdirect.co.uk](mailto:jobstdelivered@daylongdirect.co.uk) or **0808 196 6250**

\*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG  
Privacy Policy: [www.daylong.co.uk](http://www.daylong.co.uk)

**Go digital and order using JOBST Online**

Bringing you an easier way to order JOBST. Order direct with payment by account or on prescription by selecting JOBST Delivered\* as the dispenser on JOBST Online.

Scan QR code to register or visit [order.jobst.com/UK](http://order.jobst.com/UK)



\*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG  
Privacy Policy: [www.daylong.co.uk](http://www.daylong.co.uk)



JOBST®, an Essity brand



**Essity**  
T/A BSN medical Limited  
PO Box 258 • Willerby • Hull • HU10 6WT

[www.jobst.co.uk](http://www.jobst.co.uk)  
Tel: 01482 670100 • Fax: 01482 670111  
E-mail: [compression.uk@jobst.com](mailto:compression.uk@jobst.com)



By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018, the UK's implementation of the (EU) General Data Protection Regulation [GDPR]. For more information and our Privacy Policy, visit [www.jobst.co.uk](http://www.jobst.co.uk).