


# JOBST® FarrowWrap® 4000 A-D Below Knee

Ready-to-wear order form and prescription request

Email order to Customer Services on: **compression.uk@jobst.com**  
 Fax order to Customer Services on: **0345 122 3450**

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_

## Legpiece



**B-D**

Left:	Leg Length	Left:	Widest calf Circumference
Right:		Right:	

Left:	Narrowest ankle Circumference	Left:
Right:		Right:

Left:	Shoe Size	Left:	Midfoot Circumference
Right:		Right:	

Measurement only required when ordering additional JOBST FarrowHybrid (Liners)

Measurement only required when ordering additional JOBST FarrowHybrid (Liners)

## Measuring instructions

### Legpiece

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **B** (narrowest part of the ankle) to **D** (two finger widths below the crease of the knee)
3. Determine size (**XS-L**) & length (**Regular or Tall**) using chart below

When ordering JOBST FarrowWrap 4000, 1 x JOBST FarrowHybrid (Liners) wide is also provided. If additional or alternative JOBST FarrowHybrid (Liners) sizes are required, please order separately

### JOBST FarrowHybrid (Liners)

1. Measure **A1** circumference
2. Record **shoe size**
3. Determine size (**S-XL**), width (**Regular or Wide**) & colour (**Tan or Black**)

Do not apply JOBST FarrowHybrid (Liners) if calf circumference exceeds 60cm. When a large JOBST FarrowWrap 4000 is ordered, an additional non-compressive liner will be provided to accommodate a calf circumference greater than 60cm. Consider a JOBST® FarrowWrap® footpiece in this circumstance

### Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST FarrowWrap



**JOBST® FarrowWrap® 4000 Legpiece**

**Size:**  X Small  Small  Medium  Large

**Length** **B-D**  Regular (30-33 cm)  Tall (34-37 cm)

**Colour:**  Tan  Black

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

**JOBST® FarrowHybrid® (Liners)**

**Size:**  Small (2-5)  Medium (6-8)  Large (9-11)  X Large (12-13) (shoe size)

**Width** **A1**  Regular (20-27cm)  Wide (28-35cm)

**Colour:**  Tan  Black (available in Wide variant only)

PIP Code:  -

Product Code:  -

Quantity Required:

Repeat prescription required every  months

	X Small	Small	Medium	Large
<b>JOBST FarrowWrap 4000 Legpiece</b>				
<b>C</b> Calf	25-38cm	32-48cm	37-56cm	43-67cm
<b>B</b> Ankle	17-23cm	19-28cm	25-37cm	29-43cm
<b>B-D</b> Regular	30-33cm	30-33cm	30-33cm	30-33cm
<b>B-D</b> Tall	34-37cm	34-37cm	34-37cm	34-37cm

NOTE: for **B-D** length measure posteriorly (back of leg) and follow contours

# JOBST FarrowWrap 4000 Order Codes

## JOBST® FarrowWrap® 4000 Legpiece (30-40mmHg)

For patients with **mild to moderate** oedema, with or without a leg ulcer

Product Code	Size	Length	Colour	PIP Code
76661-00	Extra small	Regular	Beige	408-1543
76661-08	Extra small	Regular	Black	408-1584
76661-04	Extra small	Tall	Beige	408-1501
76661-12	Extra small	Tall	Black	408-1469
76661-01	Small	Regular	Beige	408-1550
76661-09	Small	Regular	Black	408-1592
76661-05	Small	Tall	Beige	408-1519
76661-13	Small	Tall	Black	408-1477
76661-02	Medium	Regular	Beige	408-1568
76661-10	Medium	Regular	Black	408-1600
76661-06	Medium	Tall	Beige	408-1527
76661-14	Medium	Tall	Black	408-1485
76661-03	Large	Regular	Beige	408-1576
76661-11	Large	Regular	Black	408-1618
76661-07	Large	Tall	Beige	408-1535
76661-15	Large	Tall	Black	408-1493



## JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle

Product Code	Size	Colour	PIP Code
76664-00	Small - Regular (shoe size 2-5)	Tan	408-1626
76664-01	Medium - Regular (shoe size 6-8)	Tan	408-1634
76664-02	Large - Regular (shoe size 9-11)	Tan	408-1642
76664-03	Extra Large - Regular (shoe size 12-13)	Tan	408-1659
76664-12	Small - Wide (shoe size 2-5)	Tan	408-1667
76664-13	Medium - Wide (shoe size 6-8)	Tan	408-1675
76664-14	Large - Wide (shoe size 9-11)	Tan	408-2384
76664-15	Extra Large - Wide (shoe size 12-13)	Tan	408-1683
76664-16	Small - Wide (shoe size 2-5)	Black	408-6765
76664-17	Medium - Wide (shoe size 6-8)	Black	408-6757
76664-18	Large - Wide (shoe size 9-11)	Black	408-6740
76664-19	Extra Large - Wide (shoe size 12-13)	Black	408-6732



Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm, consider a JOBST FarrowWrap footpiece in this circumstance

### Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for [codefinder.jobst.com](#) on your mobile device.



### Fast and free home delivery of prescription garments with JOBST Delivered\*

To use this dispensing service, send this form and prescription to: **JOBST Delivered, FREEPOST, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG** or order online at [order.jobst.com/UK](#) and select JOBST Delivered as the dispenser (note: prescription to be posted to JOBST Delivered unless using EPS).

If the patient is using an Electronic Prescription Service (EPS), select code **FJ708** to nominate Daylong Direct as their dispenser of choice.

For envelopes, please contact: [jobstdelivered@daylongdirect.co.uk](mailto:jobstdelivered@daylongdirect.co.uk) or **0808 196 6250**

\*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG  
Privacy Policy: [www.daylong.co.uk](#)

### Go digital and order using JOBST Online

Bringing you an easier way to order JOBST. Order direct with payment by account or on prescription by selecting JOBST Delivered\* as the dispenser on JOBST Online.

Scan QR code to register or visit [order.jobst.com/UK](#)



\*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG

Privacy Policy: [www.daylong.co.uk](#)