


JOBST® FarrowWrap®

A-G Thigh High

Ready-to-wear order form and prescription request

Email order to Customer Services on: compression.uk@jobst.com
 Fax order to Customer Services on: **0345 122 3450**

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____



Customer Service: **0345 122 3600**
 Email: compression.uk@jobst.com
 Website: www.jobst.co.uk

Measurement only required when ordering additional JOBST FarrowHybrid (Liners)

Measuring instructions

Thighpiece

1. Measure **E** (5cm above the patella) and **G** (2cm below gluteal fold) circumferences
2. Measure length medially (inside leg) following the contour from **E** to **G**
3. Determine size (**XS-XL**) & length (**Short, Regular or Tall**) using chart below
 A knee piece & a pair of appropriately sized TG Soft Liners will be provided free of charge with each JOBST FarrowWrap thighpiece

Legpiece

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **A** (floor) to **D** (two finger widths below the crease of the knee)
3. Determine size (**XS-XL**) & length (**Regular or Tall**) using chart below.
 A pair of appropriately sized liners will be provided free of charge with each JOBST FarrowWrap legpiece

Footpiece

1. Measure **A1** circumference
2. Measure length **X** from heel to base of little toe
3. Determine size (**XS-XL**) & length (**Regular or Long**) using chart below

JOBST FarrowHybrid (Liners)

1. JOBST FarrowHybrid (Liners) can be used instead of a footpiece, providing calf circumference does not exceed 60cm
2. Measure **A1** circumference
3. Record shoe size
4. Determine size (**S-XL**), width (**Regular or Wide**) & colour (**Tan or Black**)

Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST FarrowWrap



	X Small	Small	Medium	Large	X Large
Thighpiece (including kneepiece)					
G Groin	60-70cm	65-75cm	70-80cm	75-90cm	85-100cm
E Above knee	47-52cm	53-59cm	60-66cm	67-73cm	74-80cm
E-G Short	16-20cm	16-20cm	16-20cm	16-20cm	16-20cm
E-G Regular	21-26cm	21-26cm	21-26cm	21-26cm	21-26cm
E-G Tall	27-31cm	27-31cm	27-31cm	27-31cm	27-31cm

NOTE: for **E-G** length measure medially (inside leg) and follow contours

	X Small	Small	Medium	Large	X Large
Legpiece					
C Mid-Calf	36-43cm	42-50cm	48-58cm	53-63cm	58-68cm
B Ankle	21-25cm	25-30cm	30-36cm	36-42cm	42-50cm
A-D Regular	33-37cm	35-39cm	37-41cm	39-43cm	39-43cm
A-D Tall	38-41cm	40-43cm	42-45cm	44-47cm	44-47cm

NOTE: for **A-D** length measure posteriorly (back of leg) and follow contours

	X Small	Small	Medium	Large	X Large
Footpiece					
A1 Mid-foot	22-24cm	25-27cm	28-30cm	31-34cm	35-40cm
X Regular	16-17cm	17-18cm	19-20cm	20-21cm	22-23cm
X Long	18-19cm	19-20cm	21-22cm	22-23cm	24-25cm

NOTE: for **X** length measure from the heel to the base of the little toe

Note: 3cm allowance on either side of measurement

Thighpiece

Size: X Small Small Medium Large X Large

Length **E-G** Short Regular Tall **Fabric:** Lite Strong Classic

PIP Code: [][] - [][][][]

Product Code: [][][][] - [][][]

Quantity Required: []

Repeat prescription required every [] months

Legpiece

Size: X Small Small Medium Large X Large

Length **A-D** Regular Tall **Fabric:** Lite Strong Classic

PIP Code: [][] - [][][][]

Product Code: [][][][] - [][][]

Quantity Required: []

Repeat prescription required every [] months

Footpiece

Size: X Small Small Medium Large X Large

Length **X** Regular Long **Fabric:** Lite Strong Classic

PIP Code: [][] - [][][][]

Product Code: [][][][] - [][][]

Quantity Required: []

Repeat prescription required every [] months

JOBST® FarrowHybrid® (Liners)

Size: Small (2-5) Medium (6-8) Large (9-11) X Large (12-13) (shoe size)

Width **A1** Regular (20-27cm) Wide (28-35cm)

Colour: Tan Black (available in Wide variant only)

PIP Code: [][] - [][][][]

Product Code: [][][][] - [][][]

Quantity Required: []

Repeat prescription required every [] months

JOBST FarrowWrap Order Codes

JOBST® FarrowWrap® Lite (20-30mmHg)

For patients with **mild to moderate** oedema caused by venous or lymphatic conditions, where reduced compression is required

Product Code	Size	Length	Colour	PIP Code
Thighpiece (including kneepiece)				
76655-05	Extra Small	Short	Beige	360-6084
76655-00	Extra Small	Regular	Beige	360-6134
76655-10	Extra Small	Tall	Beige	360-6183
76655-06	Small	Short	Beige	360-6092
76655-01	Small	Regular	Beige	360-6142
76655-11	Small	Tall	Beige	360-6191
76655-07	Medium	Short	Beige	360-6100
76655-02	Medium	Regular	Beige	360-6159
76655-12	Medium	Tall	Beige	360-6209
76655-08	Large	Short	Beige	360-6118
76655-03	Large	Regular	Beige	360-6167
76655-13	Large	Tall	Beige	360-6217
76655-09	Extra Large	Short	Beige	360-6126
76655-04	Extra Large	Regular	Beige	360-6175
76655-14	Extra Large	Tall	Beige	360-6225
Legpiece				
76654-00	Extra Small	Regular	Beige	360-6233
76654-05	Extra Small	Tall	Beige	360-6282
76654-01	Small	Regular	Beige	360-6241
76654-06	Small	Tall	Beige	360-6290
76654-02	Medium	Regular	Beige	360-6258
76654-07	Medium	Tall	Beige	360-6308
76654-03	Large	Regular	Beige	360-6266
76654-08	Large	Tall	Beige	360-6316
76654-04	Extra Large	Regular	Beige	360-6274
76654-09	Extra Large	Tall	Beige	360-6324
Footpiece				
76652-05	Extra Small	Regular	Beige	360-6332
76652-00	Extra Small	Long	Beige	360-6381
76652-06	Small	Regular	Beige	360-6340
76652-01	Small	Long	Beige	360-6399
76652-07	Medium	Regular	Beige	360-6357
76652-02	Medium	Long	Beige	360-6407
76652-08	Large	Regular	Beige	360-6365
76652-03	Large	Long	Beige	360-6415
76652-09	Extra Large	Regular	Beige	360-6373
76652-04	Extra Large	Long	Beige	360-6423



JOBST® FarrowWrap® Strong (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for skin folds and shape distortion

Product Code	Size	Length	Colour	PIP Code
Thighpiece (including kneepiece)				
76649-05	Extra Small	Short	Beige	382-0222
76649-00	Extra Small	Regular	Beige	382-0172
76649-10	Extra Small	Tall	Beige	382-0123
76649-06	Small	Short	Beige	382-0255
76649-01	Small	Regular	Beige	382-0206
76649-11	Small	Tall	Beige	382-0156
76649-07	Medium	Short	Beige	382-0248
76649-02	Medium	Regular	Beige	382-0198
76649-12	Medium	Tall	Beige	382-0149
76649-08	Large	Short	Beige	382-0230
76649-03	Large	Regular	Beige	382-0180
76649-13	Large	Tall	Beige	382-0131
76649-09	Extra Large	Short	Beige	382-0214
76649-04	Extra Large	Regular	Beige	382-0164
76649-14	Extra Large	Tall	Beige	382-0115
Legpiece				
76648-00	Extra Small	Regular	Beige	382-0073
76648-05	Extra Small	Tall	Beige	382-0024
76648-01	Small	Regular	Beige	382-0107
76648-06	Small	Tall	Beige	382-0057
76648-02	Medium	Regular	Beige	382-0099
76648-07	Medium	Tall	Beige	382-0040
76648-03	Large	Regular	Beige	382-0081
76648-08	Large	Tall	Beige	382-0032
76648-04	Extra Large	Regular	Beige	382-0065
76648-09	Extra Large	Tall	Beige	382-0016
Footpiece				
76646-04	Extra Small	Regular	Beige	381-9976
76646-09	Extra Small	Long	Beige	381-9927
76646-05	Small	Regular	Beige	382-0008
76646-00	Small	Long	Beige	381-9950
76646-06	Medium	Regular	Beige	381-9992
76646-01	Medium	Long	Beige	381-9943
76646-07	Large	Regular	Beige	381-9984
76646-02	Large	Long	Beige	381-9935
76646-08	Extra Large	Regular	Beige	381-9968
76646-03	Extra Large	Long	Beige	381-9919



JOBST® FarrowWrap® Classic (30-40mmHg)

For patients with **moderate to severe** oedema caused by venous or lymphatic conditions, suitable for irregular shaped limbs and deep skin folds

Product Code	Size	Length	Colour	PIP Code
Thighpiece (including kneepiece)				
76658-06	Extra Small	Short	Beige	360-6431
76658-00	Extra Small	Regular	Beige	360-6480
76658-12	Extra Small	Tall	Beige	360-6530
76658-07	Small	Short	Beige	360-6449
76658-01	Small	Regular	Beige	360-6498
76658-13	Small	Tall	Beige	360-6548
76658-08	Medium	Short	Beige	360-6456
76658-02	Medium	Regular	Beige	360-6506
76658-14	Medium	Tall	Beige	360-6555
76658-09	Large	Short	Beige	360-6464
76658-03	Large	Regular	Beige	360-6514
76658-15	Large	Tall	Beige	360-6563
76658-10	Extra Large	Short	Beige	360-6472
76658-04	Extra Large	Regular	Beige	360-6522
76658-16	Extra Large	Tall	Beige	360-6571
Legpiece				
76657-00	Extra Small	Regular	Beige	360-6589
76657-05	Extra Small	Tall	Beige	360-6639
76657-01	Small	Regular	Beige	360-6597
76657-06	Small	Tall	Beige	360-6647
76657-02	Medium	Regular	Beige	360-6605
76657-07	Medium	Tall	Beige	360-6654
76657-03	Large	Regular	Beige	360-6613
76657-08	Large	Tall	Beige	360-6662
76657-04	Extra Large	Regular	Beige	360-6621
76657-09	Extra Large	Tall	Beige	360-6670
Footpiece				
76656-05	Extra Small	Regular	Beige	360-6688
76656-00	Extra Small	Long	Beige	360-6738
76656-06	Small	Regular	Beige	360-6696
76656-01	Small	Long	Beige	360-6746
76656-07	Medium	Regular	Beige	360-6704
76656-02	Medium	Long	Beige	360-6753
76656-08	Large	Regular	Beige	360-6712
76656-03	Large	Long	Beige	360-6761
76656-09	Extra Large	Regular	Beige	360-6720
76656-04	Extra Large	Long	Beige	360-6779

JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle

Product Code	Size	Colour	PIP Code
76664-00	Small - Regular (shoe size 2-5)	Tan	408-1626
76664-01	Medium - Regular (shoe size 6-8)	Tan	408-1634
76664-02	Large - Regular (shoe size 9-11)	Tan	408-1642
76664-03	Extra Large - Regular (shoe size 12-13)	Tan	408-1659
76664-12	Small - Wide (shoe size 2-5)	Tan	408-1667
76664-13	Medium - Wide (shoe size 6-8)	Tan	408-1675
76664-14	Large - Wide (shoe size 9-11)	Tan	408-2384
76664-15	Extra Large - Wide (shoe size 12-13)	Tan	408-1683
76664-16	Small - Wide (shoe size 2-5)	Black	408-6765
76664-17	Medium - Wide (shoe size 6-8)	Black	408-6757
76664-18	Large - Wide (shoe size 9-11)	Black	408-6740
76664-19	Extra Large - Wide (shoe size 12-13)	Black	408-6732

Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm, consider a JOBST FarrowWrap footpiece in this circumstance



Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for **codefinder.jobst.com** on your mobile device.



Fast and free home delivery of prescription garments with JOBST Delivered*

To use this dispensing service, send this form and prescription to: **JOBST Delivered, FREEPOST, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG** or order online at order.jobst.com/UK and select JOBST Delivered as the dispenser (note: prescription to be posted to JOBST Delivered unless using EPS).

If the patient is using an Electronic Prescription Service (EPS), select code **FJ708** to nominate Daylong Direct as their dispenser of choice.

For envelopes, please contact: jobstdelivered@daylongdirect.co.uk or **0808 196 6250**

*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG

Privacy Policy: www.daylong.co.uk

Go digital and order using JOBST Online

Bringing you an easier way to order JOBST. Order direct with payment by account or on prescription by selecting JOBST Delivered* as the dispenser on JOBST Online.

Scan QR code to register or visit order.jobst.com/UK



*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG

Privacy Policy: www.daylong.co.uk



JOBST®, an Essity brand



Essity
T/A BSN medical Limited
PO Box 258 • Willerby • Hull • HU10 6WT

www.jobst.co.uk
Tel: 01482 670100 • Fax: 01482 670111
E-mail: compression.uk@jobst.com



By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the Data Protection Act 2018, the UK's implementation of the (EU) General Data Protection Regulation [GDPR]. For more information and our Privacy Policy, visit www.jobst.co.uk.