



JOBST® Ready-to-Wear Lower Limb Compression Garments

For the management of mild to moderate venous
and lymphatic conditions



Availability on Prescription

Please use this information to check that your required JOBST garment is available on prescription before completing a prescription request form

JOBST® Opaque											
		Knee high		Thigh high						Tights	
Colour	RAL Class	Open Toe	Closed Toe	Open Toe Sensitive Top Band	Closed Toe Sensitive Top Band	Open Toe Lace Top Band	Closed Toe Lace Top Band	Open Toe Dotted Top Band	Closed Toe Dotted Top Band	Open Toe	Closed Toe
Black	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Natural	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Caramel	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Bronze	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Navy	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

JOBST® forMen Ambition		
		Knee high
Colour	RAL Class	Closed Toe
Black	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓
Dark Grey	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓
Navy	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓
Khaki	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓
Brown	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓

JOBST® UltraSheer											
		Knee high		Thigh high						Tights	
Colour	RAL Class	Open Toe	Closed Toe	Open Toe Sensitive Top Band	Closed Toe Sensitive Top Band	Open Toe Lace Top Band	Closed Toe Lace Top Band	Open Toe Dotted Top Band	Closed Toe Dotted Top Band	Open Toe	Closed Toe
Black	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	X	✓	✓	✓	X	✓	X	X	X	✓
Natural	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	X	✓	✓	✓	X	✓	X	X	X	✓
Caramel	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	X	✓	✓	✓	X	✓	X	X	X	✓
Bronze	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	X	✓	X	✓	X	✓	X	X	X	✓
Espresso	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	X	✓	X	✓	X	✓	X	X	X	✓
Anthracite	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	X	✓	X	✓	X	✓	X	X	X	✓

JOBST® forMen Explore		
		Knee high
Colour	RAL Class	Closed Toe
Black	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓
Dark Grey	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓
Navy	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓
Khaki	CCL1 (18-21mmHg) CCL2 (23-32mmHg)	✓

JOBST® Classic				
		Knee high	Thigh high	Tights
Colour	RAL Class	Open Toe	Open Toe Dotted Top Band	Open Toe
Sand	CCL2 (23-32mmHg) CCL3 (34-46mmHg)	✓	✓	✓
Pearl	CCL2 (23-32mmHg) CCL3 (34-46mmHg)	✓	✓	✓

To find out more about our extensive JOBST compression range, visit www.jobst.co.uk



JOBST® Ready-to-Wear

Prescription Request Form

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____

Measured By: _____ Tel: _____ Email: _____

Delivery Address: _____ Invoice Address: _____

Email order to Customer Services on: **compression.uk@jobst.com**

Fax order to Customer Services on: **0345 122 3450**

1. <input type="checkbox"/> JOBST® Opaque	<input type="checkbox"/> JOBST® UltraSheer	<input type="checkbox"/> JOBST® ForMen Ambition <input type="checkbox"/> JOBST® ForMen Explore	<input type="checkbox"/> JOBST® Classic
2. <input type="checkbox"/> Knee High <input type="checkbox"/> Thigh High <input type="checkbox"/> Tights	<input type="checkbox"/> Knee High <input type="checkbox"/> Thigh High <input type="checkbox"/> Tights	<input checked="" type="checkbox"/> Knee High	<input type="checkbox"/> Knee High <input type="checkbox"/> Thigh High <input type="checkbox"/> Tights
3. <input type="checkbox"/> CCL 1 (18-21mmHg) <input type="checkbox"/> CCL 2 (23-32mmHg)	<input type="checkbox"/> CCL 1 (18-21mmHg) <input type="checkbox"/> CCL 2 (23-32mmHg)	<input type="checkbox"/> CCL 1 (18-21mmHg) <input type="checkbox"/> CCL 2 (23-32mmHg)	<input type="checkbox"/> CCL 2 (23-32mmHg) <input type="checkbox"/> CCL 3 (34-46mmHg)
4. <input type="checkbox"/> Open Toe <input type="checkbox"/> Closed Toe	<input type="checkbox"/> Open Toe <input type="checkbox"/> Closed Toe	<input checked="" type="checkbox"/> Closed Toe	<input checked="" type="checkbox"/> Open Toe
5. Silicone Top Band (Thigh High only) <input type="checkbox"/> Sensitive <input type="checkbox"/> Sensitive Wide <input type="checkbox"/> Lace <input type="checkbox"/> Lace Wide <input type="checkbox"/> Dotted <input type="checkbox"/> Dotted Wide	Silicone Top Band (Thigh High only) <input type="checkbox"/> Sensitive <input type="checkbox"/> Sensitive Wide <input type="checkbox"/> Lace <input type="checkbox"/> Lace Wide <input type="checkbox"/> Dotted <input type="checkbox"/> Dotted Wide		Silicone Top Band (Thigh High only) <input type="checkbox"/> Dotted
6. <input type="checkbox"/> Petite <input type="checkbox"/> Regular	<input type="checkbox"/> Petite <input type="checkbox"/> Regular	<input type="checkbox"/> Regular <input type="checkbox"/> Long	<input type="checkbox"/> Petite <input type="checkbox"/> Regular
7. <input type="checkbox"/> Black <input type="checkbox"/> Natural <input type="checkbox"/> Caramel <input type="checkbox"/> Bronze <input type="checkbox"/> Navy	<input type="checkbox"/> Black <input type="checkbox"/> Natural <input type="checkbox"/> Caramel <input type="checkbox"/> Bronze <input type="checkbox"/> Espresso <input type="checkbox"/> Anthracite	<input type="checkbox"/> Black <input type="checkbox"/> Dark Grey <input type="checkbox"/> Navy <input type="checkbox"/> Khaki <input type="checkbox"/> Brown (Ambition only)	<input type="checkbox"/> Pearl <input type="checkbox"/> Sand


8. Size 1/I Size 2/II Size 3/III Size 4/IV Size 5/V Size 6/VI

PIP Code: <input type="text"/> <input type="text"/> <input type="text"/> - <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>	Essity Code: <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> - <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>
Quantity Required: <input type="text"/>	Repeat prescription required every <input type="text"/> months

Contact Us:
 Telephone: **0345 122 3600**
 Email: **compression.uk@jobst.com**
 Website: **www.jobst.co.uk**

Measuring Guide

JOBST® Opaque / JOBST® UltraSheer / JOBST® Classic



Measurement point Circumference (c) in cm			I	II	III	IV	V	VI	Measurements	
									Left leg	Right leg
cG (tights only)	Tights	Thigh High	43-57	45-62	49-67	53-72	56-77	60-81		
cG (silicone band wide)			49-57	53-62	57-67	61-72	65-77	69-81		
cG (silicone band)			43-48	45-52	49-56	53-60	56-64	60-68		
cF			39-52	41-56	44-60	47-65	50-69	53-73		
cE			30-37	33-40	35-43	37-45	39-48	41-51		
cD			27-33	29-36	32-39	34-42	36-45	38-48		
cC			28-34	30-37	33-40	35-43	37-46	39-49		
cBl			23-27	24-29	26-32	29-35	31-37	33-39		
cB			18-20	20-22	22-24	24-26	26-28	28-30		
cY			26-31	28-33	29-35	31-37	32-38	33-40		
cA	17-22	19-24	21-26	23-29	25-32	27-34				
			Length (I)							

ID = floor to two fingers width below crease of the knee
 IG = floor to gluteal fold
 IK = floor to crotch

Length (I) in cm:
 Knee High (ID): Petite 34-39, Regular 39-45
 Tights (IK): Petite 65-75, Regular 75-85

Thigh High (IG): Petite 60-69, Regular 70-80
 *JOBST Classic thigh high garment available in regular length only

How to use the measuring chart

Step 1

Measure the points shown on the chart and write the measurements in the left and right leg columns.

Example:

VI	Measurements	
	Left leg	Right
60-81	60	
69-81	78	


Step 2

Look along the row and circle the number ranges that your measurement fits into. This will probably be more than one column.

Example: If cA measurement is 22cm

Measurement point Circumference (c) in cm			I	II	III	IV	V	VI	M
cG (tights only)	Tights	Thigh High	43-57	45-62	49-67	53-72	56-77	60-81	
cG (silicone band wide)			49-57	53-62	57-67	61-72	65-77	69-81	
cG (silicone band)			43-48	45-52	49-56	53-60	56-64	60-68	
cF			39-52	41-56	44-60	47-65	50-69	53-73	
cE			30-37	33-40	35-43	37-45	39-48	41-51	
cD			27-33	29-36	32-39	34-42	36-45	38-48	
cC			28-34	30-37	33-40	35-43	37-46	39-49	
cBl			23-27	24-29	26-32	29-35	31-37	33-39	
cB			18-20	20-22	22-24	24-26	26-28	28-30	
cY			26-31	28-33	29-35	31-37	32-38	33-40	
cA	17-22	19-24	21-26	23-29	25-32	27-34			

JOBST® forMen Ambition / JOBST® forMen Explore



Measurement point Circumference (c) in cm			1	2	3	4	5	6	Measurements	
									Left leg	Right leg
cD	Tights	Thigh High	28-34	31-38	34-41	37-46	40-49	42-52		
cC			30-36	32-40	34-43	38-47	40-50	43-53		
cBl			23.5-27.5	26-30.5	28.5-33	31-36.5	33.5-39	35.5-41.5		
cB			18-20	20-22	22-24	24-26	26-28	28-30		
cY			26-31	28-33	29-35	31-37	32-38	33-40		
cA			17-22	19-24	21-26	23-29	25-32	27-34		
					Length (I)					

Length (I) in cm: Knee High (ID): Regular 39-45, Long 45-51

Step 3

Use the circling process for the remaining measuring points up the limb.

Example:

cA 22cm
 cY 33cm
 cB 23cm

Measurement point Circumference (c) in cm			I	II	III
cG (tights only)	Tights	Thigh High	43-57	45-62	49-67
cG (silicone band wide)			49-57	53-62	57-67
cG (silicone band)			43-48	45-52	49-56
cF			39-52	41-56	44-60
cE			30-37	33-40	35-43
cD			27-33	29-36	32-39
cC			28-34	30-37	33-40
cBl			23-27	24-29	26-32
cB			18-20	20-22	22-24
cY			26-31	28-33	29-35
cA	17-22	19-24	21-26		

Step 4

At completion, the column with the most circles is the best fit. The number at the top of this column is your size.

Example:

Measurement point Circumference (c) in cm			II	III	IV	V	VI	M
cG (tights only)	Tights	Thigh High	45-62	49-67	53-72	56-77	60-81	
cG (silicone band wide)			53-62	57-67	61-72	65-77	69-81	
cG (silicone band)			45-52	49-56	53-60	56-64	60-68	
cF			41-56	44-60	47-65	50-69	53-73	
cE			33-40	35-43	37-45	39-48	41-51	
cD			29-36	32-39	34-42	36-45	38-48	3F
cC			30-37	33-40	35-43	37-46	39-49	3F
cBl			24-29	26-32	29-35	31-37	33-39	2F
cB			20-22	22-24	24-26	26-28	28-30	2F
cY			28-33	29-35	31-37	32-38	33-40	3F
cA	19-24	21-26	23-29	25-32	27-34	2F		

JOBST® Delivered

To use this dispensing service, please send this form and prescription to:
JOBST Delivered, FREEPOST, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG.

If using an Electronic Prescription Service (EPS), please select code **FJ708** to nominate Daylong Direct as the dispenser or order online at jobstonline.co.uk and select JOBST Delivered as the dispenser.

For envelopes, please contact: jobstdelivered@daylongdirect.co.uk or **0808 196 6250**
 JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG.

By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the General Data Protection Regulation. For more information and our Privacy Policy, visit www.jobst.co.uk



JOBST®, an Essity brand



Essity
 T/A BSN medical Limited
 PO Box 258 • Willerby • Hull • HU10 6WT

www.jobst.co.uk
 Tel: 01482 670100 • Fax: 01482 670111
 E-mail: compression.uk@jobst.com



6198/122

Please note

Some patients may fit into more than one size, use clinical judgement to select most appropriate size. If the best fit is not clear, then the patient may require a made to measure flat knit garment or an adjustable compression wrap device.