

JOBST® Confidence

Custom-Fit Upper Limb Order Form

Fax order to customer services on: **0845 122 3450**
 Email order to customer services on: **compression.uk@jobst.com**

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____



JOBST Confidence

Compression Class (RAL)

CCL 1 (15-21mmHg)
 CCL 2 (23-32mmHg)

Colour

Beige Black
 Caramel Red heather
 Anthracite heather Jeans heather

Style
Armsleeve (wrist to axilla)

C-G1 Armsleeve (bias top)

Quantity

Left Right


Options

Silicone band CF
 2.5cm on top
 2.5cm inside (half silicone band)
 SoftFit
 Elbow functional zone
 Extra angled elbow (45°)*
 Deco line
 Initials (max 2 letters)

*Standard angle is 25°

Remarks

Length measurements taken on inside of arm



JOBST® Confidence Ordering Information

Compression Class / Option Code Pieces

JOBST® Confidence Armsleeve




Basic style

| | | |
|-------------------|----------|---|
| CCL 1 (15-21mmHg) | U1-16-07 | 1 |
| CCL 2 (23-32mmHg) | U2-17-07 | 1 |

Options

| | | |
|-----------------------|---------|---|
| Elbow Functional Zone | U-A008C | 1 |
| Silicone Band CF | U-A002C | 1 |
| SoftFit | U-A018C | 1 |

Options




Elbow Functional Zone

Circumferential functional zone allows greater movement and flexibility at the elbow




Extra Angled Elbow

45° angle knitted at the elbow. The standard knitted angle is 25° and is default on all armsleeves. The extra angle may be required when working continuously with a bent elbow. Standard elbow



Deco Line

A decorative seam line can be added to the front of the garment on request



Silicone Band CF

Dotted silicone band



SoftFit

Soft and secure top band that has been uniquely knitted into the garment top band for a more comfortable secure fit



Initials

Patient initials can be knitted into the garment below the JOBST logo for a personal touch