



**JOBST**

Gentle compression  
for the management of  
**mild oedema**

**JOBST® UlcerCARE™**

Ready-to-wear Compression Liners

# JOBST® UlcerCARE™

## Ready-to-Wear Compression Liners




The **innovative**, gentle answer for **mild** compression

JOBST® UlcerCARE™ Liners provide mild compression for non-ambulatory patients. Available in seven ready-to-wear sizes, they will fit an ankle circumference of up to 39cms, providing a nominal compression of 15-20mmHg\*

- ✓ Contains silk
  - to promote healthy skin
  - to facilitate donning
  - to encourage maximum concordance
- ✓ **Liner may be worn 24 hours a day**
- ✓ **Washable at 60°C** for hygienic re-use
- ✓ Helps hold wound dressing in place
- ✓ Liners available in white as standard
- ✓ Latex free

## Size Guide

Size	S	M	L	XL	2XL	3XL	4XL
Ankle	18 - 21 cm	21 - 24 cm	24 - 27 cm	27 - 30 cm	30 - 33 cm	33 - 36 cm	36 - 39 cm
Calf	29 - 35 cm	34 - 41 cm	40 - 47 cm	46 - 54 cm	48 - 56 cm	50 - 58 cm	52 - 60 cm



## Ordering Information

Order Code	Size	Pack Contents	PIP Code
73632-21	S	3	364-7666
73632-22	M	3	364-7674
73632-23	L	3	364-7682
73632-24	XL	3	364-7708
73632-25	2XL	3	364-1123
73632-26	3XL	3	364-1115
73632-27	4XL	3	364-1107

\*The nominal compression for an average ankle size



JOBST®,  
an Essity brand



4294/219

Essity  
T/A BSN medical Limited  
PO Box 258 • Willerby • Hull • HU10 6WT

www.bsnmedical.co.uk  
Tel: 01482 670100 • Fax: 01482 670111  
E-mail: compression.uk@essity.com



@JOBSTUK